

INOMax DS_{IR}[®]

(Delivery System)

Pocket Guide

Series 3 software



Automated Pre-Use Procedure
Integrated Pneumatic Backup INOMAX[®] Delivery
Transport Regulator/Cap Assembly
Oxygen Dilution Chart
INOMAX Cylinder Duration Chart
Circuit Connection Diagrams
Changing INOMAX Cylinders
High Calibration Connection Diagrams
INOMax DS_{IR} Disposable Adapters



For 24 Hour Assistance
Call 1-877-566-9466

Part No. 20751 Rev-01
2014-08

IMPORTANT: This guide is provided as a convenience and for general information only. Do not use this product without clearly and thoroughly understanding the most recent revision of the INOmax DS_{IR}[®] Operation Manual. The Operation Manual is the source for specific, updated information regarding warnings, cautions, checklists, diagrams, and/or instructions contained in this guide.

Mallinckrodt, the "M" brand mark and the Mallinckrodt
Pharmaceuticals logo are trademark of a Mallinckrodt company.
Other brands are trademarks of a Mallinckrodt company or their
respective owner. ©2016 Mallinckrodt

Contents

Automated Pre-Use Checkout.....	5
Integrated Pneumatic Backup INOMAX Delivery	8
Transport Regulator/Cap Assembly.....	10
Oxygen Dilution Chart	14
INOMAX Cylinder 88-size Duration Chart	15
INOMAX Cylinder D-size Duration Chart.....	16
Connection to Various Breathing Systems	17
Acutronics Medical Systems AG Fabian +nCPAP Evolution	18
Acutronics Medical Systems AG Fabian HFO	19
A-Plus Medical Babi-Plus Bubble CPAP.....	20
Bagging Systems While Using the Injector Module.....	22
Bunnell Life Pulse High Frequency Ventilator Circuit	24
Connecting INOmax DS _{IR} Sample Tee to the Bunnell Life Pulse Circuit	28
Connecting INOmax DS _{IR} Injector Module to the Bunnell Life Pulse Circuit	28
Circle Anesthesia System.....	29
Dräger Babylog VN500/Infinity Acute Care System and Heinen & Löwenstein Leoni-plus Ventilator	30
Fisher & Paykel Healthcare Bubble CPAP	31
Fisher & Paykel Healthcare Infant Circuit Nasal Cannula	32
Fisher & Paykel Healthcare Optiflow Breathing Circuit	33
Hamilton Arabella Nasal CPAP	34
ICU Ventilator Circuit	35
Sensormedics 3100A/B High Frequency Oscillatory Ventilator with a Filtered Circuit	36
Sensormedics 3100A/B High Frequency Oscillatory Ventilator with a Rigid or Flexible Circuit	38
SLE Life Support SLE5000.....	39
Spontaneous Breathing Patient on a Mask Circuit	40
Spontaneous Breathing Patient on a Nasal Cannula	41
Teleflex Medical Comfort Flo Humidification System.....	42
Transport Ventilator Circuit	43
Vapotherm 2000i	46
Vapotherm Precision Flow	47
Viasys Infant Flow CPAP System; Cardinal Airlife nCPAP System.....	48
Viasys Infant Flow SiPAP	49
INOblender Circuit Connection.....	50
Fisher & Paykel Healthcare Neopuff Resuscitator	51
Changing INOMAX Cylinders.....	52
High Calibration Connection Diagrams	56
INOmax DS _{IR} Disposable Adapters.....	58


Pre-Use Checkout

Automated Pre-Use Checkout

1. Turn INOmax DS_{IR} ON, verify speaker function.

Note: A low range calibration automatically starts following the self test.


A Pre-Use wizard will be displayed on the main screen, which will provide step-by-step instructions to complete the automated Pre-Use procedure.



- Pressing the NEXT button initiates the Pre-Use wizard.
- Pressing the CANCEL button exits the Pre-Use wizard. If you cancel the pre-use wizard, the manual pre-use checkout procedure can be found on the INOmax DS_{IR} Plus Pre-Use card, or in the INOmax DS_{IR} Operation Manual, section 9/Appendix.

The Pre-Use wizard can also be initiated by entering the menu screen and selecting the Pre-Use Checkout button.

Pre-Use Checkout



Pre-Use Checkout

(Intentionally left blank)

Integrated Pneumatic Backup INOMAX Delivery

Note: Use the integrated pneumatic backup function only for a short time, until a replacement delivery system can be obtained. The INOblender can also be used as a backup.

If the main delivery system fails, the INOMax DS_{IR} has an integrated pneumatic backup delivery function that allows the patient to remain connected to the ventilator. Backup NO delivery does not rely on the operation of the main system.

The INOMax DS_{IR} backup function:

- Uses a pneumatic on/off switch and a restrictor built into the delivery side of the system.
- Provides a fixed flow of INOMAX (250 mL/min) into the injector module. This fixed flow provides 20 ppm of NO when added to a continuous ventilator gas flow of 10 L/min.

WARNING:

When using the integrated pneumatic backup with breathing circuit gas flows of 5 L/min, the delivered NO dose will be approximately 40 ppm. Breathing circuit gas flows less than 5 L/min will deliver an NO dose greater than 40 ppm.

The table below indicates the nominal concentrations delivered for different ventilator gas flows.

Ventilator/Gas Flow (L/min)	5	7.5	10	15	20
NO Concentration (ppm)	40	27	20	13	10


INOMAX cylinder conc. X 0.25 L/min / ventilator flow = estimated dose

Integrated Pneumatic Backup INOMAX Delivery


The main screen:

- Indicates that backup delivery is on and that the set dose is turned off.
- Displays the estimated dose that the patient should be receiving, based on the ventilator flow.
- Displays the NO concentration table

Note: If the injector module is not functioning, the estimated backup dose graphic will be inactive.



Backup delivery mode ON (with a Low Priority alarm).



Backup delivery mode OFF.

Transport Regulator/ Cap Assembly

Transport Regulator/Cap Assembly

WARNING:

- A new INOMAX cylinder and regulator must be purged before use to ensure the patient does not receive an excess level of NO₂.
- Loss of communication between the INOmax DS_{IR} and the INOMAX cylinder for more than one hour will result in interruption of INOMAX delivery.

Caution: When using the Transport Regulator/Cap Assembly (PN 10022) ensure the cap is in place on the cylinder and the infrared cable is connected to the infrared connector port on the back of the INOmax DS_{IR}.

Note: Check the INOMAX cylinder for the correct product identity labels, cylinder concentration and expiration date. Ensure the INOMAX gas cylinder has more than 500 psig.

Step One

Note: Ensure the white plastic tip is in place on the regulator connector and not chipped or cracked (see Figure 2).

Connect a high pressure regulator to an INOMAX cylinder and tighten the fitting to the INOMAX cylinder (see Figure 1).



Figure 1




Figure 2

Transport Regulator/ Cap Assembly

Step Two

Connect the INOMAX regulator hose to one of the INOMAX inlets on the back of the INOMax DS_{IR} (see Figure 3).



Figure 3

Step Three

Connect the Infrared cable from the Transport Regulator/Cap Assembly to the back of the INOMax DS_{IR} (see Figure 4).



Figure 4

Transport Regulator/ Cap Assembly

Step Four


Place the Cap Assembly over the INOmeter (see Figure 5).

Note: Be sure to align the keyway inside the Cap Assembly with the iButton on the INOmeter (see Figure 5 and 6).



Figure 5

iButton



The electrical cord exits the cap directly above the iButton keyway

Figure 6

Step Five

Grasp the Cap Assembly to open cylinder valve (see Figure 7 and 8).



Figure 7



Figure 8

Oxygen Dilution Chart

For delivery with 800 ppm cylinder of INOMAX (nitric oxide) for inhalation.

(Illustrative Only)

		Set FiO ₂				
		.21	.40	.60	.80	1.00
INOMAX Dose	10	0.21	0.40	0.59	0.79	0.99
	20	⚠ 0.20	0.39	0.59	0.78	0.98
	40	⚠ 0.20	0.38	0.57	0.76	0.95
	80	⚠ 0.19	0.36	0.54	0.72	0.90
		Actual FiO ₂				

⚠ Caution FiO₂ less than 21%

Please Note:

The calculations on this chart have been determined based on an 800 ppm cylinder of INOMAX (nitric oxide) for inhalation.

This chart is representative of a range of doses available on the INOmax DS_{IR} and doses higher than 20 ppm are not intended as the recommended therapeutic dose.

Calculations are considered estimates and may vary under clinical conditions.


All numbers have been rounded to the nearest hundredth.

**Duration Chart
(88-size)**

INOMAX Cylinder 88-Size

For an **88-Size** 800 ppm Cylinder Concentration*
(Illustrative Only)

		FLOW			
		5 L/min	10 L/min	20 L/min	40 L/min
INOMAX Dose (ppm)	5	39 Days	19.5 Days	9.8 Days	4.9 Days
	10	19.4 Days	9.7 Days	4.8 Days	2.4 Days
	20	9.6 Days	4.8 Days	2.4 Days	1.2 Days
	40	4.7 Days	2.3 Days	1.2 Days	14 Hours
	80	2.2 Days	1.1 Days	13.3 Hours	6.6 Hours



This chart is representative of a range of doses available on the INOMax DS_{IR} and doses higher than 20 ppm are not intended as the recommended therapeutic dose.

* All calculations in the table above are based on a full cylinder, 138 bar (2000 psig), 1963 liter "88" cylinder, with a cylinder change at 14 bar (200 psig). The figures are calculated based on a total continuous breathing circuit gas flow and a cylinder conversion factor of 14.2 liters per bar/0.98 liters per psig.


- INOMAX flow = [Desired dose × total ventilator flow] ÷ [Cylinder concentration - desired dose]
- Cylinder volume = Cylinder conversion factor × cylinder pressure (bar/psig)
- Cylinder duration (hours) = (Cylinder volume ÷ INOMAX flow rate) ÷ 60

Calculations are considered estimates and may vary under clinical circumstances. For more information, call 1-877-KNOW-INO (1-877-566-9466).

**Duration Chart
(D-size)**

INOMAX Cylinder D-Size

		FLOW			
		5 L/min	10 L/min	20 L/min	40 L/min
INOMAX Dose (ppm)	5	7.0 Days	3.5 Days	1.8 Days	21 Hours
	10	3.5 Days	1.7 Days	21 Hours	10.5 Hours
	20	1.7 Days	20.7 Hours	10.3 Hours	5.2 Hours
	40	20 Hours	10 Hours	5 Hours	2.5 Hours
	80	9.5 Hours	4.8 Hours	2.4 Hours	1.2 Hours



Typically used in transport

This chart is representative of a range of doses available on the INOMax DS_{IR} and doses higher than 20 ppm are not intended as the recommended therapeutic dose.

* All calculations in the table above are based on a full cylinder, 138 bar (2000 psig), 353 liter "D" cylinder, with a cylinder change at 14 bar (200 psig). The figures are calculated based on a total continuous breathing circuit gas flow and a cylinder conversion factor of 2.6 liters per bar/0.18 liters per psig.

- INOMAX flow = [Desired dose × total ventilator flow] ÷ [Cylinder concentration - desired dose]
- Cylinder volume = Cylinder conversion factor × cylinder pressure (bar/psig)
- Cylinder duration (hours) = (Cylinder volume ÷ INOMAX flow rate) ÷ 60

Calculations are considered estimates and may vary under clinical circumstances. For more information, call 1-877-KNOW-INO (1-877-566-9466).

Circuit Connection Diagrams

Proper use of these products depends on careful reading and understanding of labeling and instructions. Please refer to the INOmax DS_{IR} and INOblender operation manuals for guidance. Also refer to the specific breathing device operation manual or instructions for use.

INOmax DS_{IR} Warnings:


- **INOmax DS_{IR} subtracts gas from the breathing circuit via the gas sampling system at 230 mL per minute; this can affect the sensitivity of a flow triggered synchronized breath mode of some ventilators. The trigger sensitivity of the ventilator should be checked after connecting the INOmax DS_{IR} to the breathing circuit.**
- **Patient disconnect and high-pressure alarms are required for the ventilator.**

INOmax DS_{IR} Cautions:

- Insert the Injector Module on the dry side of the breathing circuit prior to the humidifier (this will ensure correct flow measurement).
- Avoid medications interfering with the gas monitoring system; administer any aerosolized medications distal to the sampling tee.

**Circuit Connection
Diagrams**


**Acutronics Medical Systems AG Fabian
+nCPAP Evolution (validated for use outside
of the United States)**



- | | |
|--|--|
| 1. Fabian+ nCPAP Evolution | 8. 22F X 15M Adapter |
| 2. Patient Gas Sample Line with Nafion | 9. Humidifier |
| 3. Injector Module Electrical Cable | 10. Inspiratory Breathing Circuit Hose |
| 4. INOmax DSIR | 11. Gas Sample Tee |
| 5. NO/N ₂ Injector Tube | 12. Patient Wye |
| 6. Connecting Tube (15 inches) | 13. Proximal Pressure Tube |
| 7. Injector Module | 14. Expiratory Breathing Circuit Hose |

**Circuit Connection
Diagrams**


**Acutronics Medical Systems AG Fabian HFO
(validated for use outside of the United States)**



- | | |
|--|--|
| 1. Fabian HFO Ventilator | 9. 22F X 15M Adapter |
| 2. Patient Gas Sample Line with Nafion | 10. Humidifier |
| 3. Injector Module Electrical Cable | 11. Inspiratory Breathing Circuit Hose |
| 4. INOmax DS _{IR} | 12. Gas Sample Tee |
| 5. NO/N ₂ Injector Tube | 13. Patient Wye |
| 6. Injector Module | 14. Proximal Pressure Tube |
| 7. T-Connector Assembly, #7209.e | 15. Expiratory Breathing Circuit Hose |
| 8. Connecting Tube (15 inches) | |

**Circuit Connection
Diagrams**

A-Plus Medical Babi-Plus Bubble CPAP



- | | |
|----------------------------------|---|
| 1. Oxygen Source | 9. Tee Adapter |
| 2. Oxygen Tubing | 10. Breathing Circuit |
| 3. Pressure Relief Manifold | 11. Humidifier |
| 4. Injector Module | 12. NO/N ₂ Injector Tube |
| 5. Temperature Probe | 13. Injector Module Electrical Cable |
| 6. 90 Degree Sample Port Adapter | 14. INOmax DS _{IR} |
| 7. Nasal Prongs | 15. Patient Gas Sample Line with Nafion |
| 8. Babi Plus Bubble PAP Valve | |

**Circuit Connection
Diagrams**


(Intentionally left blank)

Bagging Systems While Using the Injector Module

WARNING: To minimize the delivered concentration of NO₂, the following steps should be taken for use with the manual resuscitator bags:

- Use the smallest bag adequate to deliver the desired tidal volume.
- Oxygen tubing lengths greater than 72 inches should not be used.
- Use the highest fresh gas flow rate (up to 15 L/min) that is practical.
- Use the lowest practical inspired oxygen concentration.
- After starting fresh gas flow, squeeze the bag several times to empty residual gas in the bag prior to using the system to ventilate a patient.

Circuit Connection Diagrams



1. O₂ Flowmeter (wall outlet or cylinder)
2. O₂ Tubing
3. 15M X 4.5 mm Adapter
4. 22M/15F X 22M/15F Adapter
5. Injector Module
6. 15M X 4.5 mm Adapter
7. O₂ Tubing
8. O₂ Tubing Sample Tee
9. Patient Gas Sample Line with Nafion
10. NO/N₂ Injector Tube
11. Resuscitator Bag with O₂ Reservoir
12. Injector Module Electrical Cable


Circuit Connection Diagrams

WARNING:

- The hyperinflation bag will, under some conditions, contain NO₂ in excess of 1 ppm. Use of large tidal volume breaths may expose the patients to the NO₂ present in the bag, for part of the breath. In general, if the inspiratory flow rate induced by the manual ventilation does not exceed the fresh gas flow rate, the patient should not be exposed to the concentrations of NO₂ present in the hyperinflation bag.
- Adult and infant hyperinflation bags generate more NO₂ when used at lower minute ventilation. If use of the bag is interrupted (for example, to adjust a tracheal tube), before resuming ventilation of the patient, the user should squeeze the bag several times to empty residual gas from the bag.
- Because of the potential for inhalation of excessive concentrations of NO₂, and the difficulty in monitoring the peak inhaled NO₂ concentrations, ventilation with a hyperinflation bag or self inflating bag is intended only for short term use.
- The monitoring system within the INOmax DSIR will not detect generation of NO₂ within the hyperinflation bag or self-inflating bag devices and the alarms for excessive NO₂ cannot warn of NO₂ produced within the manual bag system.
- To minimize the delivered concentration of NO₂, the following steps should be taken for use with the manual resuscitator bags:
 - Concentrations greater than 20 ppm NO should not be used because of excessive NO₂ generation.
 - Use the smallest bag adequate to deliver the desired tidal volume.
 - Oxygen tubing lengths greater than 72 inches should not be used.
 - Use the highest fresh gas flow rate (up to 15 L/min) that is practical.
 - Use the lowest practical inspired oxygen concentration.
 - After starting fresh gas flow, squeeze the bag several times to empty residual gas in the bag prior to using the system to ventilate a patient.

**Circuit Connection
Diagrams**

**Bagging Systems While Using the Injector
Module**



- | | |
|--|-------------------------------|
| 1. O ₂ Flowmeter | 8. Pressure Gauge |
| 2. Injector Module Electrical Cable | 9. 15M X 4.5mm Adapter |
| 3. NO/N ₂ Injector Tube | 10. Injector Module |
| 4. Patient Gas Sample Line with Nafion | 11. 22M/15F X 22M/15F Adapter |
| 5. O ₂ Tubing | 12. 15M X 4.5mm Adapter |
| 6. O ₂ Tubing Sample Tee | 13. O ₂ Tubing |
| 7. Hyper-Inflation Bag | |

Circuit Connection Diagrams

Bunnell Life Pulse High Frequency Ventilator Circuit

WARNING:


- The INOmax DS_{IR} backup mode (250 mL/min.) should not be used with the Bunnell Life Pulse as ventilator flow rates are normally below the recommended ventilator flows.
- Place the Life Pulse in Standby prior to suctioning the patient to avoid NO delivery transiently exceeding the set dose by up to 30 ppm. Press ENTER to reestablish ventilation as soon as the catheter is removed from the airway. This will limit the extent of over delivery above the NO set dose.

Caution:

- If the set dose is below 5 ppm and the Servo pressure is 2.0 psig. or less, this will result in flow rates outside of the specification of the Injector Module and fluctuating NO values may result.
- A one-way valve should be placed between the injector module and the humidifier chamber to prevent water from backing up into the injector module if the Life Pulse is either put into Standby or cycled OFF.
- There are higher pressures in the breathing circuit than normal; use only parts provided in disposable package #50248 and tightly secure all connections.

Circuit Connection
Diagrams

**Bunnell Life Pulse
High Frequency Ventilator Circuit (cont.)**




- | | |
|----------------------------|---|
| 1. INOmax DS _{IR} | 8. Sample Tee |
| 2. Bunnell Life Pulse | 9. Patient Box |
| 3. Humidifier | 10. One-Way Valve |
| 4. Humidifier | 11. Injector Module |
| 5. Conventional Ventilator | 12. NO/N ₂ Injector Tube |
| 6. Life Port Adapter | 13. Injector Module Electrical Cable |
| 7. Endotracheal Tube | 14. Patient Gas Sample Line with Nafion |


**Circuit Connection
Diagrams**

**Connecting INOmax DS_{IR} Sample Tee to the
Bunnell Life Pulse Circuit**

1. From Patient Box
2. Cut Green tube at midpoint (approximately six inches from the Life Port Adapter)
3. Patient Gas Sample Line with Nafion
4. Insert Sample Tee
5. Life Port Adapter
6. Endotracheal Tube




**Connecting INOmax DS_{IR} Injector Module to the
Bunnell Life Pulse Circuit**



- | | |
|------------------------------|-------------------------------------|
| 1. Gas Out Tube from Vent | 5. 15M X 4.5mm I.D. Adapter |
| 2. 15M X 4.5mm I.D. Adapter | 6. 3cm Piece of Green Gas Out Tube |
| 3. 22M/15F X 22M/15F Adapter | 7. One-Way Valve |
| 4. Injector Module | 8. Green Gas Out Tube to Humidifier |

Circuit Connection Diagrams

Circle Anesthesia System




- | | |
|---|-------------------------------|
| 1. Patient Gas Sample Line with Nafion | 8. Absorber Inspiratory Port |
| 2. Patient Gas Sample Line Input Connection | 9. Absorber |
| 3. INOmax DS _{IR} | 10. Injector Module |
| 4. Bellows Assembly | a. Injector Module Input End |
| 5. Ventilator | b. Injector Module Output End |
| 6. Ventilator Drive Gas | 11. Inspiratory Tubing |
| 7. Absorber Expiratory Port | 12. 22M/15F X 22M/15F Adapter |
| | 13. Gas Sample Tee |
| | 14. Patient Wye |

WARNING: Fresh gas flow should be equal to or greater than patient minute ventilation to avoid recirculation of gases.

**Circuit Connection
Diagrams**


**Dräger Babylog VN500/Infinity Acute Care
System and Heinen & Löwenstein
Leoni-plus Ventilator (validated for use outside
of the United States)**



- | | |
|---|--|
| 1. Patient wye | 9. Injector Module |
| 2. Dräger Babylog VN500 / Leoni-plus Ventilator | 10. One-Way Valve |
| 3. Ventilator Expiratory Port | 11. Humidifier Inlet |
| 4. Ventilator Inspiratory Port | 12. Humidifier |
| 5. Patient Gas Sample Line Input Connection | 13. Humidifier Outlet |
| 6. INOmax DS _{IR} | 14. Patient Gas Sample Line
with Nafion |
| 7. NO/N ₂ Injector Tube Front Panel Connection | 15. Gas Sample Tee |
| 8. Injector Module Electrical Cable Front Panel
Connection | |

Circuit Connection
Diagrams


Fisher & Paykel Bubble CPAP



- | | |
|----------------------------------|---|
| 1. Oxygen Source | 10. F/P Inline Infant Nebulizer Kit (RT010) Adapter |
| 2. Oxygen Tubing | 11. Breathing Circuit |
| 3. Bubble CPAP Pressure Manifold | 12. Humidifier |
| 4. 22F X 15M Adapter | 13. NO/N ₂ Injector Tube |
| 5. 22M/15F X 22M/15F Adapter | 14. Injector Module Electrical Cable |
| 6. Injector Module | 15. INOmax DS _{1R} |
| 7. Temperature Probe | 16. Patient Gas Sample Line with Nafion |
| 8. Nasal Prong Infant Interface | |
| 9. Bubble CPAP Generator | |

**Circuit Connection
Diagrams**


Fisher & Paykel Infant Circuit Nasal Cannula



- | | |
|--|-------------------------------------|
| 1. Patient Gas Sample Line with Nafion | 9. Injector Module Electrical Cable |
| 2. INOmax DS _{IR} | 10. NO/N ₂ Injector Tube |
| 3. Oxygen Source | 11. Humidifier |
| 4. Oxygen Tubing | 12. Breathing Circuit |
| 5. 22F X 15M Adapter | 13. Temperature Probe |
| 6. Injector Module | 14. Gas Sample Tee |
| 7. Pressure Relief Manifold | 15. Nasal Cannula |
| 8. 22M/15F X 22M/15F Adapter | |

**Circuit Connection
Diagrams**


Fisher & Paykel Optiflow Breathing Circuit



- | | |
|--|--------------------------------------|
| 1. Patient Gas Sample Line with Nafion | 10. Breathing Circuit |
| 2. INOmax DS _{IR} | 11. Temperature Probe |
| 3. Oxygen Source | 12. Gas Sample Tee |
| 4. Breathing Circuit Hose | 13. 22M/15F X 22M/15F Adapter |
| 5. Injector Module | 14. 22 mm ID X 22 mm ID Cuff Adapter |
| 6. Injector Module Electrical Cable | 15. Optiflow Tracheostomy |
| 7. NO/N ₂ Injector Tube | 16. Optiflow Nasal Cannula |
| 8. 22F X 15M Adapter | 17. Optiflow Mask |
| 9. Humidifier | |

**Circuit Connection
Diagrams**


Hamilton Arabella Nasal CPAP



- | | |
|--|-----------------------------------|
| 1. Arabella | 8. Humidifier |
| 2. Patient Gas Sample Line with Nafion | 9. Heated Delivery Circuit |
| 3. INOmax DSIR | 10. Temperature Probe |
| 4. NO/N ₂ Injector Tube | 11. Universal Generator |
| 5. Injector Module Electrical Cable | 12. Arabella Sample Tee |
| 6. Injector Module | 13. 90 Degree Sample Port Adapter |
| 7. 22F X 15M Adapter | |


**Circuit Connection
Diagrams**

ICU Ventilator Circuit



**Circuit Connection
Diagrams**

**Sensormedics 3100A/B High Frequency
Oscillatory Ventilator with a Filtered Circuit**




Circuit Connection Diagrams

- | | |
|---|---|
| 1. Sensormedics 3100A/B Ventilator | 10. Paw Limit Valve Control |
| 2. Ventilator Outlet | 11. Filter |
| 3. 22M Adapter | 12. Humidifier Inlet |
| 4. Injector Module | 13. Humidifier Outlet |
| 5. Injector Module Electrical Cable
Connection | 14. Bias Flow Tube |
| 6. INOmax DS _{IR} | 15. Patient Gas Sample Line with Nafion |
| 7. NO/N ₂ Injector Tube | 16. 90 Degree Sample Port Adapter |
| 8. 8 mm Tubing X 15M Adapter | 17. Dump Valve Control |
| 9. One-Way Valve | 18. Paw Control Valve |

WARNING: Omission of the one-way valve may result in high NO delivery.

**Circuit Connection
Diagrams**

**Sensormedics 3100A/B High Frequency Oscillatory
Ventilator with a Rigid or Flexible Circuit**



- | | |
|--|---|
| 1. Sensormedics 3100A/B Ventilator | 9. Humidifier Inlet |
| 2. Ventilator Outlet | 10. Humidifier Outlet |
| 3. 22M Adapter | 11. Patient Gas Sample Line with Nafion |
| 4. Injector Module | 12. 90 Degree Sample Port Adapter |
| 5. INOmax DS _{IR} | 13. Bias Flow Tube |
| 6. NO/N ₂ Injector Tube Connection | |
| 7. Injector Module Electrical Cable Connection | |
| 8. One-Way Valve | |


WARNING: Omission of the one-way valve may result in high NO delivery.

Circuit Connection Diagrams

SLE Life Support SLE5000

Note: • Validated for use outside of the United States.


- A one-way valve is not required for use during high frequency ventilation mode.



- | | |
|--|--|
| 1. Patient Wye | 7. Injector Module |
| 2. Patient Gas Sample Line with Nafion | 8. 22F X 15M Adapter |
| 3. SLE5000 | 9. Humidifier |
| 4. INOmax DS _{IR} | 10. Inspiratory Breathing Circuit Hose |
| 5. NO/N ₂ Injector Tube | 11. Expiratory Breathing Circuit Hose |
| 6. Injector Module Electrical Cable | 12. Gas Sample Tee |

**Circuit Connection
Diagrams**

**Spontaneous Breathing Patient on a Mask
Circuit**




- | | |
|------------------------------|--|
| 1. O ₂ Tubing | 10. One-Way Valve |
| 2. 15M X 4.5 mm Adapter | 11. Sealed Face Mask |
| 3. 22M/15F X 22M/15F Adapter | 12. One-Way Valve |
| 4. Breathing Circuit Tee | 13. Patient Gas Sample Line with Nafion |
| 5. Breathing Circuit Bag | 14. NO/N ₂ Injector Tube |
| 6. Injector Module | 15. INOmax DSIR |
| 7. Breathing Circuit Hose | 16. Injector Module Electrical Cable |
| 8. Gas Sample Tee | 17. O ₂ Flowmeter (wall outlet or cylinder) |
| 9. 22M/15F X 22M/15F Adapter | |

Circuit Connection Diagrams

Spontaneous Breathing Patient on a Nasal Cannula

The INOmax DS_{IR} can be used with nasal cannula to deliver INOMAX concentrations from 5-80 ppm and an oxygen flow rate as low as 2 L/min.


WARNING: Do not use the INOmax DS_{IR} backup mode with flow rates less than 5 L/min.



- | | |
|------------------------------|--|
| 1. O ₂ Flowmeter | 7. O ₂ Tubing Sample Tee |
| 2. O ₂ Tubing | 8. Patient Nasal Cannula |
| 3. 15M x 4.5 mm Adapter | 9. Patient Gas Sample Line with Nafion |
| 4. 22M/15F x 22M/15F Adapter | 10. INOmax DS _{IR} |
| 5. 300 mm of 22 mm Hose | 11. Injector Module Electrical Cable |
| 6. Injector Module | 12. NO/N ₂ Injector Tube |

**Circuit Connection
Diagrams**


Teleflex Medical Comfort Flo Humidification System



- | | |
|---|--|
| 1. Patient Gas Sample Line with Nafion | 10. Temperature Probe Connector |
| 2. INOmax DS _{IR} | 11. Second Temperature Probe Connector |
| 3. Injector Module | 12. Comfort Flo Cannula |
| 4. System Pressure Relief Valve | 13. Injector Module Electrical Cable |
| 5. Air/Oxygen Blender or Oxygen Blender | 14. NO/N ₂ Injector Tube |
| 6. Oxygen Tubing | 15. 22F X 15M Adapter |
| 7. Temperature Probe (Short Cable) | 16. ConchaTherm Heated Humidifier |
| 8. Angled 22 mm Connector | 17. Temperature Probe (Long Cable) |
| 9. Patient Circuit | 18. 90 Degree Sample Port Adapter |

**Circuit Connection
Diagrams**


Transport Ventilator Diagram



- | | |
|--|--|
| 1. Patient Wye | 7. Ventilator Inspiratory Port |
| 2. Expiratory Breathing Circuit Hose | 8. 22M/15F X 22M/15F Adapter |
| 3. Patient Gas Sample Line with Nafion | 9. Injector Module Electrical Cable |
| 4. Ventilator Expiratory Valve | 10. NO/N ₂ Injector Tube |
| 5. Ventilator | 11. Injector Module |
| 6. INOmax DS _{IR} | 12. Inspiratory Breathing Circuit Hose |
| | 13. Gas Sample Tee |

**Circuit Connection
Diagrams**

Single-Limb Transport Ventilator Diagram




- | | |
|--|--|
| 1. PEEP/Exhalation Valve | 7. Ventilator Inspiratory Port |
| 2. Patient Wye | 8. 22M/15F X 22M/15F Adapter |
| 3. Circuit Hose | 9. Injector Module Electrical Cable |
| 4. Patient Gas Sample Linewith
Nafion | 10. NO/N ₂ Injector Tube |
| 5. Ventilator | 11. Injector Module |
| 6. INOmax DS _{IR} | 12. Inspiratory Breathing Circuit Hose |
| | 13. Gas Sample Tee |

**Circuit Connection
Diagrams**

(Intentionally left blank)

**Circuit Connection
Diagrams**

Vapotherm 2000i




- | | |
|------------------------------|--|
| 1. INOmax DS _{IR} | 9. 15M x 4.5mm Adapter |
| 2. O ₂ Flowmeter | 10. Vapotherm 2000i |
| 3. O ₂ Tubing | 11. Patient Delivery Tube |
| 4. 15M x 4.5mm Adapter | 12. O ₂ Tubing Sample Tee |
| 5. 22M/15F x 22M/15F Adapter | 13. Patient Cannula |
| 6. 300mm of 22mm Hose | 14. Patient Gas Sample Line
with Nafion |
| 7. 22M/15F x 22M/15F Adapter | 15. NO/N ₂ Injector Tube |
| 8. Injector Module | 16. Injector Module Electrical Cable |

Circuit Connection Diagrams


Connection to the Vapotherm Precision Flow

- The INOmax DS_{IR} adds NO/N₂ gas flow to the breathing circuit flow in proportion to the NO setting (up to 10% at 80 ppm) and subtracts gas from the breathing circuit via gas sampling at a nominal flow rate of 0.23 L/min.
- These effects impact the delivered gas flow rate when using the Vapotherm Precision Flow. It is recommended that after an NO setting change the user checks the delivered gas flow rate and adjusts the gas source flow rate as necessary.
- Follow all manufacturer instructions for connection to the Vapotherm Precision Flow.



**Circuit Connection
Diagrams**


**Viasys Infant Flow CPAP System;
Cardinal AirLife nCPAP System**



- | | |
|----------------------------|--|
| 1. INOmax DS _{IR} | 7. Patient Gas Sample Line with Nafion |
| 2. Heated Delivery Circuit | 8. Humidifier |
| 3. Infant Flow System | 9. 22F X 15M Adapter |
| 4. Infant Flow Generator | 10. Injector Module |
| 5. Sample Tee | 11. NO/N ₂ Injector Tube |
| 6. Temperature Probe | 12. Injector Module Electrical Cable |

Circuit Connection
Diagrams

Viasys Infant Flow SiPAP



- | | |
|---------------------------------|---|
| 1. INOmax DS _{IR} | 8. Heated Delivery Circuit |
| 2. Abdominal Respiratory Sensor | 9. Humidifier |
| 3. Transducer Interface | 10. 22F X 15M Adapter |
| 4. Infant Flow SiPAP | 11. Injector Module |
| 5. Infant Flow Generator | 12. Injector Module Electrical Cable |
| 6. Sample Tee | 13. NO/N ₂ Injector Tube |
| 7. Temperature Probe | 14. Patient Gas Sample Line with Nafion |

INOblender Circuit Connection Diagram

INOblender Warnings:

- **The purge procedure must be followed to help ensure NO₂ is purged from the pressure regulator, INOblender and hoses before the manual resuscitator bag is connected to the patient. The manual bag should be squeezed continuously during use to avoid NO₂ building up in the bag. If the bag is not squeezed continuously while delivering INOMAX, the bag should be removed from the patient and the purge procedure performed before continuing.**
- **Persons using this device should be trained on and experienced in the use of this device to assure effective administration of INOMAX and to avoid injury to the patient or others resulting from inhalation of excess INOMAX, nitrogen dioxide or other reaction products.**


INOblender Cautions:

- Refer to the manufacturer's procedures for using the resuscitation bag. When finished, turn the INOMAX cylinder off and continue to flow O₂ until the NO pressure gauge reads zero, then turn the O₂ flow off and the NO dial to zero ppm.

Note: Connections to various ventilators as well as their corresponding disposable circuits, are unique to each manufacturer. Please refer to the specific breathing device operation manual or instructions for use for guidance.

**INOblender Circuit
Connection Diagram**

**INOblender Connection to the Fisher & Paykel
Neopuff Resuscitator**



1. Oxygen Source
2. Neopuff
3. T-Piece Circuit (with Duckbill Port)
4. Patient Connection
5. Temperature Probe
6. Humidified Resuscitation System Circuit
7. Humidifier
8. Oxygen Tubing
9. INOblender
10. INOMAX Inlet

Changing INOMAX Cylinders

Changing INOMAX Cylinders

WARNING:

- A new INOMAX cylinder and regulator must be purged before use to ensure the patient does not receive an excess level of NO₂.
- Loss of communication between the INOmax DS_{IR} and the INOMAX cylinder for more than one hour will result in interruption of INOMAX delivery.

Caution:

- Replace an INOMAX cylinder when its pressure is less than 200 psig.
- When using the Transport Regulator/Cap Assembly (PN 10022) ensure the cap is fully seated and in place on the INOmeter and the infrared cable is connected and latched to the infrared connector port on the back of the INOmax DS_{IR}.


Note: Ensure the white plastic tip is in place.

- A. Attach a regulator to an INOMAX cylinder with greater than 500 psig.**




Changing INOMAX Cylinders


**B. Perform high
pressure leak test.**



**C. Purge the high
pressure hose.**



**D. Connect the
pressure hose.**



Changing INOMAX Cylinders

Changing INOMAX Cylinders (cont.)

- E. Open the cylinder valve (this may activate the “Two Cylinders Open” alarm until the empty cylinder valve is closed).**

Note: If using the INOmax DS_{IR} Transport Regulator/Cap Assembly, transfer the cap from the exhausted INOMAX cylinder to the new INOMAX cylinder at this time; the “Cylinder Not Detected” alarm may occur.

- F. Close the cylinder valve on the empty cylinder and remove the hose from the back of the INOmax DS_{IR}.**




- G. Depressurize and remove the regulator from the empty cylinder.**

**Changing INOMAX
Cylinders**

(Intentionally left blank)

High Calibration Connection Diagrams

Connection Diagram for NO and NO₂ High Range Calibration




1. Cylinder Cap
2. Regulator Seal
3. Regulator
4. Tubing Adapter
5. Attach Tube Kit
6. Turn valve counter-clockwise to start gas flow

- 6a. If the pressure is in the red or black zone (0-25 psig) select another INOcal cylinder.
7. Attach tube kit to sample line
8. Calibrate sensor
9. Remove sample line from tube kit
10. Turn valve clockwise to stop gas flow
11. One-way Valve

High Calibration
Connection Diagrams

Calibration Setup for O₂ High Range Calibration



1. 100% O₂ Source
2. O₂ Tubing
3. 15M x 4.5mm I.D. Adapter
4. Gas Sample Tee
5. Patient Gas Sample Line with Nafion

**INomax DS_{IR}
Disposable Adapters**


INomax DS_{IR} Patient Circuit Disposables

(Note: Graphics not actual size)


Adapter,
15M Fits 4.5mm ID Tubing



Adapter, 22F X 15M




Adapter, 22M/15F X 22M/15F



Adapter, Cuff, 22mm ID X 22mm ID




Adapter, Gas Sample Tee



Adapter, 90 degree Sample Port




Bunnell Life Pulse Disposable
Adapters Convenience Pack




Disk Filter, 0.5 micron



Neonatal Tubing,
10mm (2 pieces)




NO/N₂ Injector Tube




**INOmax DS_{IR}
Disposable Adapters**

One-way Valve, 22F X 22M



Patient Gas Sample Line with Nafion



Pediatric Extension,
15 mm (6 inches)




Sample Tee, O₂ Tubing



Sensormedics 3100A/B
Filtered Circuit
Disposable Adapters
Convenience Pack



Water Separator Cartridge



**Mallinckrodt Manufacturing LLC
6603 Femrite Drive,
Madison, WI 53718-6801 USA
1-877-566-9466**

Part No. 20751 Rev-01
2014-08