



Technical Bulletin

Ref: **TB- 15001**

Product: INOmax DS_{IR}[®] / INOblender[®]

Affected parts: INOmax DS_{IR} / INOblender

Subject: **Device Label – UDI Bar Code Addition**

Issue date: **September 2015**

Priority: Low

Classification: Information for Distribution

Due to new standards required by the FDA, a Unique Device Identifier (UDI) barcode will be added to INOmax DS_{IR} and INOblender device labels. The addition of this bar code does not require the user to change their equipment management processes; however, the look of the labels will change slightly.

The components of the labels are as follows:

REF: Stock number

S/N: Serial number

UDI Components:

(01) Device identifier

(21) Serial number

The table below displays examples (final size and detail may vary slightly) of the current device labels and the new UDI device labels (which include the UDI bar code):

Device	Current Device Label	New UDI Device Label
INOmax DS _{IR}		
INOblender		

If you have further questions, please contact technical support at 1-877-566-9466 (North America) or your local representative.

Mallinckrodt Manufacturing LLC
6603 Femrite Drive,
Madison, WI 53718-6801 USA

Mallinckrodt, the "M" brand mark and the Mallinckrodt Pharmaceuticals logo are trademark of a Mallinckrodt company. Other brands are trademarks of a Mallinckrodt company or their respective owner. ©2016 Mallinckrodt